

**Error message when opening Tools 1.2 : "This form has been created to use illuminant #85712 but there are only three active illuminants. Will adjust to current illuminant!".**

Replace the topbar.cnf and botbar.cnf files, by reinstalling Datacolor Tools. These files are located in the following folders:

c:\documents and settings\all users\application data\Datacolor\Tools\Tools global

and also in sub-folders of:

c:\documents and settings\all users\application data\Datacolor

and On Vista or Windows 7:

C:\ProgramData\Datacolor\Tools\Tools\_Global

C:\ProgramData\Datacolor\Tools\Tools\_AppUser\_DCI

C:\ProgramData\Datacolor\Tools\Tools\_AppUser\_user

^

^

*Author: Datacolor*

*Saved From: <http://knowledgebase.datacolor.com/article-537.html>*